Imię

Data

Nazwisko

Kierunek studiów

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1. Axis Lock

Click to lock or unlock axis scaling.

2. Measurement Selector

Click to change the displayed measurement or units.

3. Axis

Click and drag to scale each axis individually.

4. Plot Area

Click and drag to move the plot area.

5. Legend

Select a run to analyze by clicking the run in the legend. Right-click the run symbol to change the run color and other data set properties.

6. Tools Drawer

Click to show or hide the Graph tools.

7. Scale to Fit

Automatically scale the axes to fit all data from the selected run within the plot area.

8. Cursor Mode

Click to toggle the function of the cursor.

- Switch to **Move Mode** to use the cursor to move the plot area.
- Switch to **Select Mode** to highlight an area on the plot for data analysis.

9. SPARKvue provides two types of coordinates tools:

- **Coordinates Tool**
Use to determine the coordinates of a single data point. Use the delta tool to determine the difference between two points of data.
- **Multi-Coordinates Tool**
Use to compare coordinates across multiple runs or measurements.

10. Slope Tool

The slope tool provides the slope of a line tangent to a curve at a single point.

11. Statistics

11. Statistics

Use the statistics tool to find the minimum, maximum, mean, and standard deviation of a range of data.

12. SPARKvue provides two curve fit tools:

- **Linear Fit**
Use the linear fit tool to fit a line to data and determine the slope and y-intercept of the trend line.
- **Curve Fit Tool**
Use the curve fit tool to fit a curve to data and obtain values about the trend line parameters.

13. Annotation Tool

Use the annotation tool to label a data point with a brief note.

14. Use these tools to display multiple y-axes:

- **Add y-axis**
Display multiple measurements on the same plot area by adding additional y-axes to a line graph.
- **Align Origins**
When displaying more than one y-axis, select to align the origins and scales of each axis.
- **Add Plot Area**
Display an additional measurement in a line graph on separate plot area located below the current plot area.

15. Prediction Tool


Before collecting data, draw a prediction of the results directly on the graph.

16. Properties

Change the visibility of data points, symbols, and lines connecting data points in a line graph.

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1. Axis Lock

Click to lock or unlock axis scaling.

2. Measurement Selector

Click to change the displayed measurement or units.

3. Axis

Click and drag to scale each axis individually.

4. Plot Area

Click and drag to move the plot area.

5. Legend

Select a run to analyze by clicking the run in the legend. Right-click the run symbol to change the run color and other data set properties.

6. Scale to Fit Tool

Change the scale of a line graph to zoom in on a range of data or to view an entire run of data.

7. Coordinates Tool

Use to determine the coordinates of a single data point. Use the delta tool to determine the difference between two points of data.

7. Coordinates Tool

Use to determine the coordinates of a single data point. Use the delta tool to determine the difference between two points of data.

8. T Annotation Tool

Use the annotation tool to label a data point in a graph with a brief note.

9. Normalize Data

Select to change the vertical scale to make 1 the maximum value.

10. Change Sample Rate

Click these tools to increase or decrease the sample rate.

Tip

Set a sample rate that is approximately twice the frequency range displayed on the FFT.

11. Properties

Change the visibility of data points, symbols, and lines connecting data points in a line graph.

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